

c-Abl(Phospho-Tyr412) Antibody

Catalog No: #11125

Package Size: #11125-1 50ul #11125-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


| | |
|--------------------|--------------------------------|
| Product Name | c-Abl(Phospho-Tyr412) Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Applications | IHC IF |
| Species Reactivity | Hu Ms Rt |
| Immunogen Type | Peptide-KLH |
| Target Name | c-Abl |
| Modification | Phospho |
| Other Names | ABL, JTK7, p150, c-ABL, |
| Accession No. | Swiss-Prot:P00519 Gene ID:25 |
| Uniprot | P00519 |
| GeneID | 25; |
| Calculated MW | 122kd |
| Concentration | 1.0mg/ml |

Application Details


Immunohistochemistry: 1:50~1:100

Immunofluorescence: 1:100~1:200

Images



Immunofluorescence staining of methanol-fixed HeLa cells using c-Abl (phospho-Tyr412) antibody (#11125, Red).



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using c-Abl (phospho-Tyr412) antibody (#11125).

Note: This product is for in vitro research use only