

c-Abl(Phospho-Tyr412) Antibody

Catalog No: #11125

Package Size: #11125-1 50ul #11125-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	c-Abl(Phospho-Tyr412) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	IHC IF
Species Reactivity	Hu Ms Rt
Immunogen Type	Peptide-KLH
Target Name	c-Abl
Modification	Phospho
Other Names	ABL, JTK7, p150, c-ABL,
Accession No.	Swiss-Prot:P00519 Gene ID:25
Uniprot	P00519
GeneID	25;
Calculated MW	122kd
Concentration	1.0mg/ml

Application Details


Immunohistochemistry: 1:50~1:100

Immunofluorescence: 1:100~1:200

Images



Immunofluorescence staining of methanol-fixed HeLa cells using c-Abl (phospho-Tyr412) antibody (#11125, Red).



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using c-Abl (phospho-Tyr412) antibody (#11125).

Note: This product is for in vitro research use only