

Met (Phospho-Tyr1234/1235) Antibody

Catalog No: #11585

Package Size: #11585-1 50ul #11585-2 100ul


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Met (Phospho-Tyr1234/1235) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB
Species Reactivity	Hu Ms Rt
Immunogen Type	Peptide-KLH
Target Name	Met
Modification	Phospho
Other Names	HGFR, AUTS9, RCCP2, c-Met ,
Accession No.	Swiss-Prot:Q13085Gene ID:31
Uniprot	Q13085
GeneID	31;
Concentration	1.0mg/ml

Images



Western blot analysis of extracts from HUVEC cells untreated or treated with EGF using Met (Phospho-Tyr1234/1235) Antibody #11585.

Note: This product is for in vitro research use only