

Puma Antibody

Catalog No: #21688


Package Size: #21688-1 50ul #21688-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Puma Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB
Species Reactivity	Hu
Immunogen Type	Peptide-KLH
Target Name	Puma
Other Names	BBC3, JFY1, JFY-1,
Accession No.	Swiss-Prot:Q9BXH1Gene ID:27113
Uniprot	Q96PG8
GeneID	27113;
Concentration	1.0mg/ml

Images



Western blot analysis of extracts from 293 and MCF-7 cells using Puma Antibody #21688.

Note: This product is for in vitro research use only