

GTSE-1 Antibody FITC Conjugated

Catalog No: #C04158F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GTSE-1 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 710-741 741 derived from human GTSE-1
Conjugates	FITC
Target Name	GTSE-1
Other Names	B99; G2 and S phase-expressed protein 1; GTSE-1; Protein B99 homolog; GTSE1
Accession No.	Swiss-Prot#Q9NYZ3NCBI Gene ID51512
Uniprot	Q9NYZ3
GeneID	51512;
Excitation Emission	494nm 518nm
Cell Localization	Cytoplasm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

May be involved in p53-induced cell cycle arrest in G2 M phase by interfering with microtubule rearrangements that are required to enter mitosis. Overexpression delays G2 M phase progression.

Note: This product is for in vitro research use only