

CD163 M130 Antibody FITC Conjugated

Catalog No: #C04167F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD163 M130 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 800-850 1156 derived from human CD163 M130
Conjugates	FITC
Target Name	CD163 M130
Other Names	M13; MM13; Scavenger receptor cysteine-rich type 1 protein M13; Hemoglobin scavenger receptor; CD163
Accession No.	Swiss-Prot#Q86VB7NCBI Gene ID9332
Uniprot	Q86VB7
GeneID	9332;
Excitation Emission	494nm 518nm
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1ug/Test IF=1:50-200

Background

Acute phase-regulated receptor involved in clearance and endocytosis of hemoglobin haptoglobin complexes by macrophages and may thereby protect tissues from free hemoglobin-mediated oxidative damage. May play a role in the uptake and recycling of iron, via endocytosis of hemoglobin haptoglobin and subsequent breakdown of heme. Binds hemoglobin haptoglobin complexes in a calcium-dependent and pH-dependent manner. Exhibits a higher affinity for complexes of hemoglobin and multimeric haptoglobin of HP*1F phenotype than for complexes of hemoglobin and dimeric haptoglobin of HP*1S phenotype. Induces a cascade of intracellular signals that involves tyrosine kinase-dependent calcium mobilization, inositol triphosphate production and secretion of IL6 and CSF1. Isoform 3 exhibits the higher capacity for ligand endocytosis and the more pronounced surface expression when expressed in cells. After shedding, the soluble form (sCD163) may play an anti-inflammatory role, and may be a valuable diagnostic parameter for monitoring macrophage activation in inflammatory conditions.

Note: This product is for in vitro research use only