

Versican Antibody FITC Conjugated

Catalog No: #C04170F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Versican Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 3360-3395 3396 derived from human Versican core protein
Conjugates	FITC
Target Name	Versican
Other Names	WGN; ERVR; GHAP; PG-M; WGN1; CSPG2; Versican core protein; Chondroitin sulfate proteoglycan core protein 2; Chondroitin sulfate proteoglycan 2; Glial hyaluronate-binding protein; Large fibroblast proteoglycan; VCAN
Accession No.	Swiss-Prot#P13611NCBI Gene ID1462
Uniprot	P13611
GeneID	1462;
Excitation Emission	494nm 518nm
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

May play a role in intercellular signaling and in connecting cells with the extracellular matrix. May take part in the regulation of cell motility, growth and differentiation. Binds hyaluronic acid.

Note: This product is for in vitro research use only