

IL-3 Antibody FITC Conjugated

Catalog No: #C04209F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IL-3 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	MsB RtB
Immunogen Description	KLH conjugated synthetic peptide aa 100-150 166 derived from mouse IL-3
Conjugates	FITC
Target Name	IL-3
Other Names	BPA; PSF; HCGF; Il-3; MCGF; Csfmu; Interleukin-3; Hematopoietic growth factor; Mast cell growth factor; Multipotential colony-stimulating factor; P-cell-stimulating factor; Il3
Accession No.	Swiss-Prot#P01586NCBI Gene ID16187
Uniprot	P01586
GeneID	16187;
Excitation Emission	494nm 518nm
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200B

Background

Granulocyte macrophage colony-stimulating factors are cytokines that act in hematopoiesis by controlling the production, differentiation, and function of 2 related white cell populations of the blood, the granulocytes and the monocytes-macrophages. This CSF induces granulocytes, macrophages, mast cells, stem cells, erythroid cells, eosinophils and megakaryocytes.

Note: This product is for in vitro research use only