

E2-Tag Mouse Monoclonal Antibody

Catalog No: #T609

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	E2-Tag Mouse Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Purification	Affinity-chromatography
Applications	WB
Specificity	E2- Tag antibody can recognize C-terminal, internal, and N-terminal E2-tagged proteins.
Immunogen Type	Peptide-KLH
Immunogen Description	E2-tag peptide GVSSTSSDFRDR conjugated to KLH.
Target Name	E2-Tag
SDS-PAGE MW	N/A
Concentration	1.0mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20 °C.Stable for 12 months from date of receipt

Application Details

WB: 1:3000~1:10000

Images

2ug E2 fusion protein+ Primary antibody dilution at
1 γ • 1:5,000 2 γ • 1:10,000

Note: This product is for in vitro research use only