

## IL17RA CD217 Antibody FITC Conjugated

Catalog No: #C04216F

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	IL17RA CD217 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	HuB MsB RtB
Immunogen Description	KLH conjugated synthetic peptide aa 240-280 866 derived from human IL17R
Conjugates	FITC
Target Name	IL17RA CD217
Other Names	CD217; IL17R; CANDF5; CDw217; IL-17RA; hIL-17R; Interleukin-17 receptor A; IL-17 receptor A; IL17RA
Accession No.	Swiss-Prot#Q96F46NCBI Gene ID23765
Uniprot	Q96F46
GeneID	23765;
Excitation Emission	494nm 518nm
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200B

## Background

Receptor for IL17A, IL17F and, in dimer with IL17RE, for IL17C. Binds its IL17A ligand with low affinity, suggesting that additional components are involved in IL17A-induced signaling.

Note: This product is for in vitro research use only