

IL-17RB Antibody FITC Conjugated

Catalog No: #C04219F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IL-17RB Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt
Species Reactivity	HuB MsB RtB B
Immunogen Description	KLH conjugated synthetic peptide aa 350-400 502 derived from human IL-17RB
Conjugates	FITC
Target Name	IL-17RB
Other Names	CRL4; EVI27; IL17BR; IL17RH1; Interleukin-17 receptor B; IL-17 receptor B; IL-17RB; Cytokine receptor-like 4; IL-17 receptor homolog 1; IL-17Rh1; Interleukin-17B receptor; IL-17B receptor; IL17RB; UNQ251 PRO19612
Accession No.	Swiss-Prot#Q9NRM6NCBI Gene ID55540
Uniprot	Q9NRM6
GeneID	55540;
Excitation Emission	494nm 518nm
Cell Localization	Cytoplasm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200B

Background

Receptor for the proinflammatory cytokines IL17B and IL17E. May play a role in controlling the growth and or differentiation of hematopoietic cells.

Note: This product is for in vitro research use only