

## IL-22 Antibody FITC Conjugated

Catalog No: #C04228F

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	IL-22 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt IF
Species Reactivity	HuB MsB RtB B
Immunogen Description	KLH conjugated synthetic peptide aa 33-83 179 derived from human IL-22
Conjugates	FITC
Target Name	IL-22
Other Names	TIFa; IL-21; IL-22; ILTIF; IL-TIF; IL-D110; zcyto18; TIFIL-23; Interleukin-22; Cytokine Zcyto18; IL-10-related T-cell-derived-inducible factor; IL22; UNQ3099 PRO10096
Accession No.	Swiss-Prot#Q9GZX6NCBI Gene ID50616
Uniprot	Q9GZX6
GeneID	50616;
Excitation Emission	494nm 518nm
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200B IF=1:50-200B

## Background

Cytokine that contributes to the inflammatory response in vivo.

Note: This product is for in vitro research use only