

Integrin alpha 6 Antibody HRP Conjugated

Catalog No: #C04235H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Integrin alpha 6 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WBB B IHC-P
Species Reactivity	HuB MsB RtB B B
Immunogen Description	KLH conjugated synthetic peptide aa 1060-1110 1130 derived from mouse Integrin alpha 6
Conjugates	HRP
Target Name	Integrin alpha 6
Other Names	CD49f; VLA-6; ITGA6B; Integrin alpha-6; CD49 antigen-like family member F; ITGA6
Accession No.	Swiss-Prot#P23229NCBI Gene ID3655
Uniprot	P23229
GeneID	3655;
Excitation Emission	N A
Cell Localization	Cytoplasm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000B B IHC-P=1:50-200B

Background

Integrin alpha-6 beta-1 is a receptor for laminin on platelets. Integrin alpha-6 beta-4 is a receptor for laminin in epithelial cells and it plays a critical structural role in the hemidesmosome.

Note: This product is for in vitro research use only