

CD85c Antibody FITC Conjugated

Catalog No: #C04249F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD85c Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt IF
Species Reactivity	HuB MsB Rt
Immunogen Description	KLH conjugated synthetic peptide aa 150-200 590 derived from human LIR8 CD85c
Conjugates	FITC
Target Name	CD85c
Other Names	LIR8; CD85C; LIR-8; Leukocyte immunoglobulin-like receptor subfamily B member 5; CD85 antigen-like family member C; Leukocyte immunoglobulin-like receptor 8; LILRB5
Accession No.	Swiss-Prot#O75023NCBI Gene ID10990
Uniprot	O75023
GeneID	10990;
Excitation Emission	494nm 518nm
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200B IF=1:50-200B

Background

May act as receptor for class I MHC antigens.

Note: This product is for in vitro research use only