

C5b-9 Antibody FITC Conjugated

Catalog No: #C04256F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C5b-9 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt IF
Species Reactivity	HuB MsB RtB
Immunogen Description	Purified Human SC5b-9 Complex
Conjugates	FITC
Target Name	C5b-9
Other Names	C9D; ARMD15; Complement component C9; C9
Accession No.	Swiss-Prot#P02748NCBI Gene ID735
Uniprot	P02748
GeneID	735;
Excitation Emission	494nm 518nm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=2ug/TestB IF=1:50-200

Background

Constituent of the membrane attack complex (MAC) that plays a key role in the innate and adaptive immune response by forming pores in the plasma membrane of target cells. C9 is the pore-forming subunit of the MAC.

Note: This product is for in vitro research use only