

TSHB Mouse Polyclonal Antibody HRP Conjugated

Catalog No: #C04259H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TSHB Mouse Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WBB B IHC-PB IHC-F
Species Reactivity	HuB MsB RtB B B B B
Immunogen Description	KLH conjugated synthetic peptide derived from human TSHB
Conjugates	HRP
Target Name	TSHB
Other Names	CHNG4; Thyroid stimulating hormone beta subunit; Thyroid stimulating hormone, beta precursor; Thyrotropin beta subunit; TSHB.
Accession No.	NCBI Gene ID7252
Uniprot	P01222
GeneID	7252;
Excitation Emission	N A
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000B B IHC-P=1:50-200B IHC-F=1:50-200B

Background

Thyroid stimulating hormone, also known as thyrotropin, is secreted from cells in the anterior pituitary called thyrotrophs, finds its receptors on epithelial cells in the thyroid gland, and stimulates that gland to synthesize and release thyroid hormones. TSH is a glycoprotein hormone composed of two subunits which are non covalently bound to one another. The alpha subunit of TSH is also present in two other pituitary glycoprotein hormones, follicle stimulating hormone and luteinizing hormone, and, in primates, in the placental hormone chorionic gonadotropin. Each of these hormones also has a unique beta subunit, which provides receptor specificity. In other words, TSH is composed of alpha subunit bound to the TSH beta subunit, and TSH associates only with its own receptor. Free alpha and beta subunits have essentially no biological activity.

Note: This product is for in vitro research use only