

## KLF2 Antibody FITC Conjugated

Catalog No: #C04312F

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	KLF2 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt IF
Species Reactivity	HuB MsB RtB B B B B
Immunogen Description	KLH conjugated synthetic peptide aa 265-315 355 derived from human KLF2
Conjugates	FITC
Target Name	KLF2
Other Names	LKLF; Krueppel-like factor 2; Lung krueppel-like factor; KLF2
Accession No.	Swiss-Prot#Q9Y5W3NCBI Gene ID10365
Uniprot	Q9Y5W3
GeneID	10365;
Excitation Emission	494nm 518nm
Cell Localization	Nucleus
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200B IF=1:50-200B

## Background

KLF2 Activates transcription by binding to the CACCC box in the beta-globin gene promoter.

Note: This product is for in vitro research use only