

## TNFRSF13B Antibody HRP Conjugated

Catalog No: #C04313H

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	TNFRSF13B Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WBB B IHC-PB IHC-F
Species Reactivity	MsB RtB
Immunogen Description	KLH conjugated synthetic peptide aa 20-70 249 derived from mouse TNFRSF13B
Conjugates	HRP
Target Name	TNFRSF13B
Other Names	Taci; 1200009E08Rik; Tumor necrosis factor receptor superfamily member 13B; Transmembrane activator and CAML interactor; CD267; Tnfrsf13b
Accession No.	Swiss-Prot#Q9ET35NCBI Gene ID57916
Uniprot	Q9ET35
GeneID	57916;
Excitation Emission	N A
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:500-2000B B IHC-P=1:50-200B IHC-F=1:50-200B

## Background

Receptor for TNFSF13 APRIL and TNFSF13B TALL1 BAFF BLYS that binds both ligands with similar high affinity. Mediates calcineurin-dependent activation of NF-AT, as well as activation of NF-kappa-B and AP-1. Involved in the stimulation of B- and T-cell function and the regulation of humoral immunity (By similarity).

Note: This product is for in vitro research use only