

FUCA1 Antibody FITC Conjugated

Catalog No: #C04378F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FUCA1 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	HuB MsB RtB B
Immunogen Description	KLH conjugated synthetic peptide aa 350-400 466 derived from human FUCA1
Conjugates	FITC
Target Name	FUCA1
Other Names	FUCA; Tissue alpha-L-fucosidase; Alpha-L-fucosidase I; Alpha-L-fucoside fucohydrolase 1; Alpha-L-fucosidase 1; FUCA1; Nbla10230
Accession No.	Swiss-Prot#P04066NCBI Gene ID2517
Uniprot	P04066
GeneID	2517;
Excitation Emission	494nm 518nm
Cell Localization	Extracellular, Lysosome
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200B

Background

Alpha-L-fucosidase is responsible for hydrolyzing the alpha-1,6-linked fucose joined to the reducing-end N-acetylglucosamine of the carbohydrate moieties of glycoproteins.

Note: This product is for in vitro research use only