

Cyclin E1 Antibody

Catalog No: #29030

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cyclin E1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	Cyclin E1 Antibody detects endogenous levels of total Cyclin E1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human Cyclin E1
Target Name	Cyclin E1
Other Names	alpha-fetoprotein, FETA, HPAFP
Accession No.	Swiss-Prot:P12004Gene ID:5111
Uniprot	P12004
GeneID	5111;
SDS-PAGE MW	48KD
Concentration	1.0mg/ml
Formulation	Supplied at 2mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.3, 0.05% sodium azide and 50% glycerol.
Storage	Store at -20°C/1 year

Application Details

WB: 1:500-1:2000

Background

A member of the highly conserved cyclin family, whose members are characterized by a dramatic periodicity in protein abundance through the cell cycle. Cyclins function as regulators of CDK kinases.

Note: This product is for in vitro research use only