

## B-Raf (Phospho-Ser365) Antibody HRP Conjugated

Catalog No: #C04427H

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	B-Raf (Phospho-Ser365) Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WBB B IHC-PB IHC-F
Species Reactivity	HuB MsB RtB B B
Immunogen Description	KLH conjugated synthetic phosphopeptide aa 350-390 766 derived from human B Raf around the phosphorylation site of Ser365
Conjugates	HRP
Target Name	B-Raf Ser365
Other Names	NS7; BRAF1; RAFB1; B-RAF1; Serine threonine-protein kinase B-raf; Proto-oncogene B-Raf; p94; v-Raf murine sarcoma viral oncogene homolog B1; BRAF
Accession No.	Swiss-Prot#P15056NCBI Gene ID673
Uniprot	P15056
GeneID	673;
Excitation Emission	N A
Cell Localization	Cytoplasm, Nucleus, Cell membrane
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:500-2000B B IHC-P=1:50-200B IHC-F=1:50-200B

## Background

Protein kinase involved in the transduction of mitogenic signals from the cell membrane to the nucleus. May play a role in the postsynaptic responses of hippocampal neuron. Phosphorylates MAP2K1, and thereby contributes to the MAP kinase signal transduction pathway.

Note: This product is for in vitro research use only