

CEBP alpha (Phospho-Thr226) Antibody HRP Conjugated

Catalog No: #C04456H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CEBP alpha (Phospho-Thr226) Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WBB B IHC-PB IHC-F
Species Reactivity	HuB MsB RtB B B B
Immunogen Description	KLH conjugated synthetic phosphopeptide aa 200-250 358 derived from human CEBP alpha around the phosphorylation site of Thr226
Conjugates	HRP
Target Name	CEBP alpha Thr226
Other Names	CEBP; C EBP-alpha; CCAAT enhancer-binding protein alpha; C EBP alpha; CEBPA
Accession No.	Swiss-Prot#P49715NCBI Gene ID1050
Uniprot	P49715
GeneID	1050;
Excitation Emission	N A
Cell Localization	Nucleus
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000B B IHC-P=1:50-200B IHC-F=1:50-200B

Background

C EBP is a DNA-binding protein that recognizes two different motifs: the CCAAT homology common to many promoters and the enhanced core homology common to many enhancers.

Note: This product is for in vitro research use only