

## CD64 IGFR1 Antibody PE Conjugated

Catalog No: #C04753P

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CD64 IGFR1 Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WBB
Species Reactivity	HuB MsB RtB B B
Immunogen Description	KLH conjugated synthetic peptide aa 100-150 374 derived from human CD64 IGFR1
Conjugates	PE
Target Name	CD64 IGFR1
Other Names	CD64; FCRI; CD64A; IGFR1; High affinity immunoglobulin gamma Fc receptor I; IgG Fc receptor I; Fc-gamma RI; Fc-gamma RIA; FcgammaRIa; FCGR1A; FCG1; FCGR1
Accession No.	Swiss-Prot#P12314NCBI Gene ID2209
Uniprot	P12314
GeneID	2209;
Excitation Emission	480,565nm 578nm
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:500-2000B B

## Background

High affinity receptor for the Fc region of immunoglobulins gamma. Functions in both innate and adaptive immune responses.

Note: This product is for in vitro research use only