

## CD21 EBV receptor Antibody PE Conjugated

Catalog No: #C04926P

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CD21 EBV receptor Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 1050-1092 1092 derived from human CD21
Conjugates	PE
Target Name	CD21 EBV receptor
Other Names	CR; C3DR; CD21; CVID7; SLEB9; Complement receptor type 2; Cr2; Complement C3d receptor; Epstein-Barr virus receptor; EBV receptor
Accession No.	Swiss-Prot#P20023NCBI Gene ID1380
Uniprot	P20023
GeneID	1380;
Excitation Emission	480,565nm 578nm
Cell Localization	Cytoplasm, Cell membrane
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200 IF=1:50-200

## Background

Receptor for complement C3Dd, for the Epstein-Barr virus on human B-cells and T-cells and for HNRPU. Participates in B lymphocytes activation.

Note: This product is for in vitro research use only