

## HNF4A Antibody FITC Conjugated

Catalog No: #C04956F

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |   |
|-----------------------|---|
| Product Name          | HNF4A Antibody FITC Conjugated  |
| Host Species          | Rabbit  |
| Clonality             | Polyclonal  |
| Isotype               | IgG   |
| Purification          | Purified by Protein A.  |
| Applications          | Flow-Cyt IF   |
| Species Reactivity    | Hu Ms Rt  |
| Immunogen Description | KLH conjugated synthetic peptide aa 300-330 464 derived from human HNF4 alpha   |
| Conjugates            | FITC  |
| Target Name           | HNF4A   |
| Other Names           | TCF; HNF4; MODY; FRTS4; MODY1; NR2A1; TCF14; HNF4a7; HNF4a8; HNF4a9; NR2A21; HNF4alpha; Hepatocyte nuclear factor 4-alpha; HNF-4-alpha; Nuclear receptor subfamily 2 group A member 1; Transcription factor 14; TCF-14; Transcription factor HNF-4; HNF4A |
| Accession No.         | Swiss-Prot#P41235NCBI Gene ID3172   |
| Uniprot               | P41235  |
| GeneID                | 3172;   |
| Excitation Emission   | 494nm 518nm   |
| Cell Localization     | Nucleus   |
| Concentration         | 1mg ml  |
| Formulation           | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.  |
| Storage               | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.   |

## Application Details

Flow-Cyt=1:50-200 IF=1:50-200

## Background

Transcriptionally controlled transcription factor. Binds to DNA sites required for the transcription of alpha 1-antitrypsin, apolipoprotein CIII, transthyretin genes and HNF1-alpha. May be essential for development of the liver, kidney and intestine.

Note: This product is for in vitro research use only