

Megalyn Antibody FITC Conjugated

Catalog No: #C05010F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|--|
| Product Name | Megalyn Antibody FITC Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Purification | Purified by Protein A. |
| Applications | ICC IF |
| Species Reactivity | Hu Ms Rt |
| Immunogen Description | KLH conjugated synthetic peptide aa 3411-3450 4655 derived from human Megalyn |
| Conjugates | FITC |
| Target Name | Lrp2 Megalyn |
| Other Names | DBS; GP33; Low-density lipoprotein receptor-related protein 2; LRP-2; Glycoprotein 33; Megalyn; LRP2 |
| Accession No. | Swiss-Prot#P98164NCBI Gene ID4036 |
| Uniprot | P98164 |
| GeneID | 4036; |
| Excitation Emission | 494nm 518nm |
| Cell Localization | Extracellular |
| Concentration | 1mg ml |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

ICC=1:50-200 IF=1:50-200

Background

Acts together with cubilin to mediate HDL endocytosis (By similarity). May participate in regulation of parathyroid-hormone and para-thyroid-hormone-related protein release.

Note: This product is for in vitro research use only