

ASGPR1 Antibody HRP Conjugated

Catalog No: #C05158H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ASGPR1 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 250-290 291 derived from human ASGPR1
Conjugates	HRP
Target Name	ASGPR1
Other Names	HL-1; ASGPR; ASGPR1; CLEC4H1; Asialoglycoprotein receptor 1; ASGP-R 1; ASGPR 1; C-type lectin domain family 4 member H1; Hepatic lectin H1; ASGR1
Accession No.	Swiss-Prot#P07306NCBI Gene ID432
Uniprot	P07306
GeneID	432;
Excitation Emission	N A
Cell Localization	Extracellular
Concentration	1mg/ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000

Background

Mediates the endocytosis of plasma glycoproteins to which the terminal sialic acid residue on their complex carbohydrate moieties has been removed. The receptor recognizes terminal galactose and N-acetylgalactosamine units. After ligand binding to the receptor, the resulting complex is internalized and transported to a sorting organelle, where receptor and ligand are disassociated. The receptor then returns to the cell membrane surface.

Note: This product is for in vitro research use only