

Pancreatic Lipase Antibody HRP Conjugated

Catalog No: #C05163H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pancreatic Lipase Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human Pancreatic Lipase
Conjugates	HRP
Target Name	Pancreatic Lipase
Other Names	PL; Pancreatic lipase; Pancreatic triacylglycerol lipase; pancrelipase; PTL; LIPP_HUMAN.
Accession No.	NCBI Gene ID5406
Uniprot	P16233
GeneID	5406;
Excitation Emission	N A
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Background

Pancreatic lipase, is involved in the hydrolysis of dietary triglycerides to fatty acids in the intestine. This hydrolysis, which is initiated in the stomach by gastric lipase, is essential for the intestinal absorption of long-chain triglyceride fatty acids.

Note: This product is for in vitro research use only