

MAP4K1 Antibody FITC Conjugated

Catalog No: #C05169F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MAP4K1 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human MAP4K1
Conjugates	FITC
Target Name	MAP4K1
Other Names	Hematopoietic progenitor kinase 1; Hematopoietic progenitor kinase; HPK1; Human Hematopoietic Progenitor Kinase; M4K1_HUMAN; Map4k1; MAPK ERK kinase kinase 1; MEK kinase kinase 1; MEKKK 1; Mitogen Activated Protein Kinase Kinase Kinase Kinase 1; Mitogen-activated protein kinase kinase kinase
Accession No.	NCBI Gene ID11184
Uniprot	Q92918
GeneID	11184;
Excitation Emission	494nm 518nm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Background

May play a role in the response to environmental stress. Appears to act upstream of the JUN N-terminal pathway. May play a role in hematopoietic lineage decisions and growth regulation.

Note: This product is for in vitro research use only