

BNIP3 Antibody FITC Conjugated

Catalog No: #C05242F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	BNIP3 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human BNIP3
Conjugates	FITC
Target Name	BNIP3
Other Names	BCL2 Adenovirus E1B 19kDa Interacting Protein 3; BCL2 adenovirus E1B 19 kDa protein interacting protein 3; BCL2 adenovirus E1B 19 kDa protein-interacting protein 3; BNIP 3; BNIP3; BNIP3_HUMAN; NIP 3; NIP3.
Accession No.	NCBI Gene ID664
Uniprot	Q12983
GeneID	664;
Excitation Emission	494nm 518nm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

Apoptosis-inducing protein that, which can overcome BCL2 suppression. May play a role in repartitioning calcium between the two major intracellular calcium stores in association with BCL2.

Note: This product is for in vitro research use only