

## testis antigen 1 Antibody HRP Conjugated

Catalog No: #C05290H

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	testis antigen 1 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F
Species Reactivity	Hu
Immunogen Description	KLH conjugated synthetic peptide aa 80-130 180 derived from human Cancer testis antigen 1B 1A
Conjugates	HRP
Target Name	testis antigen 1
Other Names	CTAG; ESO1; CT6.1; CTAG1; LAGE-2; LAGE2B; NY-ESO-1; Cancer testis antigen 1; Autoimmunogenic cancer testis antigen NY-ESO-1; Cancer testis antigen 6.1; L antigen family member 2; CTAG1A; LAGE2; LAGE2A
Accession No.	Swiss-Prot#P78358NCBI Gene ID1485
Uniprot	P78358
GeneID	1485;246100;
Excitation Emission	N A
Cell Localization	Cytoplasm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

## Background

The protein encoded by this gene is an antigen that is overexpressed in many cancers but that is also expressed in normal testis. This gene is found in a duplicated region of the X-chromosome and therefore has a neighboring gene of identical sequence. [provided by RefSeq, Jan 2012]

Note: This product is for in vitro research use only