

PRRSV-GP5 Antibody HRP Conjugated

Catalog No: #C05292H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PRRSV-GP5 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P
Species Reactivity	Virus
Immunogen Description	KLH conjugated synthetic peptide derived from PRRSV GP5: 101-200/200
Conjugates	HRP
Target Name	PRRSV-GP5
Other Names	Porcine reproductive and respiratory syndrome, PRRS.
Accession No.	Swiss-Prot#O55492
Uniprot	O55492
Excitation Emission	N A
Concentration	1 mg/ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Western blotting: 1:100-1000

Immunohistochemistry1:100-500

Background

GP5 glycosylated envelope protein

Note: This product is for in vitro research use only