

CD8B Antibody PE Conjugated

Catalog No: #C05446P

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD8B Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	ICC IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 80-130 210 derived from human CD8B
Conjugates	PE
Target Name	CD8B
Other Names	LY3; P37; LEU2; LYT3; CD8B1; T-cell surface glycoprotein CD8 beta chain; CD8b
Accession No.	Swiss-Prot#P10966NCBI Gene ID926
Uniprot	P10966
GeneID	926;
Excitation Emission	480,565nm 578nm
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Background

Identifies cytotoxic suppressor T-cells that interact with MHC class I bearing targets. CD8 is thought to play a role in the process of T-cell mediated killing.

Note: This product is for in vitro research use only