

human secretory IgA Antibody HRP Conjugated

Catalog No: #C05461H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	human secretory IgA Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IHC-P IHC-F
Species Reactivity	Hu
Immunogen Description	Purified human sIgA
Conjugates	HRP
Target Name	human secretory IgA
Other Names	GIF; GIFB; GRIF; ZnMT3; Metallothionein-3; MT-3; Growth inhibitory factor; Metallothionein-III; MT-III; MT3
Accession No.	Swiss-Prot#P25713NCBI Gene ID4504
Uniprot	P25713
GeneID	4504;
Excitation Emission	N A
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200

Background

Binds heavy metals. Contains three zinc and three copper atoms per polypeptide chain and only a negligible amount of cadmium. Inhibits survival and neurite formation of cortical neurons in vitro.

Note: This product is for in vitro research use only