

## nNOS(Ab-852) Antibody

Catalog No: #21159

Package Size: #21159-1 50ul #21159-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	nNOS(Ab-852) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Applications	WB IF
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous level of total nNOS protein.
Immunogen Type	Peptide-KLH
Immunogen Description	Peptide sequence around aa.850~854 (F-N-S-V-S) derived from Human nNOS.
Target Name	nNOS
Other Names	BNOS; Constitutive NOS; N-NOS; NC-NOS; NOS1
Accession No.	Swiss-Prot: P29475NCBI Protein: NP_000611.1
Uniprot	P29475
GeneID	4842;
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


## Application Details

Predicted MW: 160kd


Western blotting: 1:500~1:1000

Immunofluorescence: 1:100~1:200

## Images



Western blot analysis of extracts from Rat and Mouse brain tissue using nNOS(Ab-852) Antibody #21159.



Immunofluorescence staining of methanol-fixed HeLa cells using nNOS(Ab-852) Antibody #21159.

## Background

Produces nitric oxide (NO) which is a messenger molecule with diverse functions throughout the body. In the brain and peripheral nervous system, NO displays many properties of a neurotransmitter.

Liu JL et al. (2008) World J Gastroenterol. 14(20):3242-8.

Guggilam A et al. (2008) Eur J Heart Fail. 10(7):625-34

Marquina M et al. (2008) Br J Dermatol.159(1):68-76.

Note: This product is for in vitro research use only