

Axin 2 Antibody FITC Conjugated

Catalog No: #C05957F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|--|
| Product Name | Axin 2 Antibody FITC Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Purification | Purified by Protein A. |
| Applications | IF |
| Species Reactivity | Hu Ms Rt |
| Immunogen Description | KLH conjugated synthetic peptide aa 350-400 843 derived from human Axin 2 |
| Conjugates | FITC |
| Target Name | Axin 2 |
| Other Names | AXIL; ODCRCS; Axin-2; Axin-like protein; Axis inhibition protein 2; Conductin; AXIN2 |
| Accession No. | Swiss-Prot#Q9Y2T1NCBI Gene ID8313 |
| Uniprot | Q9Y2T1 |
| GeneID | 8313; |
| Excitation Emission | 494nm 518nm |
| Cell Localization | Cytoplasm |
| Concentration | 1mg ml |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

IF=1:50-200

Background

Inhibitor of the Wnt signaling pathway. Down-regulates beta-catenin. Probably facilitate the phosphorylation of beta-catenin and APC by GSK3B (By similarity).

Note: This product is for in vitro research use only