

## LGALS3BP Antibody HRP Conjugated

Catalog No: #C05965H

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	LGALS3BP Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 340-390 585 derived from human LGALS3BP
Conjugates	HRP
Target Name	LGALS3BP
Other Names	9K; M2BP; gp9; CyCAP; BTBD17B; MAC-2-BP; TANGO1B; Galectin-3-binding protein; Basement membrane autoantigen p15; Lectin galactoside-binding soluble 3-binding protein; Mac-2-binding protein; MAC2BP; Mac-2 BP; Tumor-associated antigen 9K; LGALS3BP
Accession No.	Swiss-Prot#Q08380NCBI Gene ID3959
Uniprot	Q08380
GeneID	3959;
Excitation Emission	N A
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

## Background

Promotes integrin-mediated cell adhesion. May stimulate host defense against viruses and tumor cells.

Note: This product is for in vitro research use only