

## LANCL2 Antibody HRP Conjugated

Catalog No: #C05988H

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	LANCL2 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 365-415 450 derived from human LANCL2
Conjugates	HRP
Target Name	LANCL2
Other Names	TASP; GPR69B; LanC-like protein 2; Testis-specific adriamycin sensitivity protein; LANCL2
Accession No.	Swiss-Prot#Q9NS86NCBI Gene ID55915
Uniprot	Q9NS86
GeneID	55915;
Excitation Emission	N A
Cell Localization	Cytoplasm, Nucleus, Cell membrane
Concentration	1mg ml
Formulation	10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% Gentamicin.
Storage	Store at 4C for 12 months.

## Application Details

Western blotting: 1:100-1000Immunohistochemistry1:100-500

## Background

Necessary for abscisic acid (ABA) binding on the cell membrane and activation of the ABA signaling pathway in granulocytes.

Note: This product is for in vitro research use only