

## H Cadherin Antibody HRP Conjugated

Catalog No: #C06029H

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	H Cadherin Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 340-390 716 derived from human H Cadherin
Conjugates	HRP
Target Name	H Cadherin
Other Names	CDHH; P105; Cadherin-13; Heart cadherin; H-cadherin; Truncated cadherin; T-cad; T-cadherin; CDH13
Accession No.	Swiss-Prot#P55290NCBI Gene ID1012
Uniprot	P55290
GeneID	1012;
Excitation Emission	N A
Cell Localization	Cell membrane
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IHC-P=1:50-200 IHC-F=1:50-200

## Background

Cadherins are calcium-dependent cell adhesion proteins. They preferentially interact with themselves in a homophilic manner in connecting cells; cadherins may thus contribute to the sorting of heterogeneous cell types. May act as a negative regulator of neural cell growth.

Note: This product is for in vitro research use only