

ADAM23 Antibody FITC Conjugated

Catalog No: #C06054F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ADAM23 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 580-630 832 derived from human ADAM23
Conjugates	FITC
Target Name	ADAM23
Other Names	MDC3; MDC-3; Disintegrin and metalloproteinase domain-containing protein 23; ADAM 23; Metalloproteinase-like, disintegrin-like, and cysteine-rich protein 3; ADAM23
Accession No.	Swiss-Prot#O75077NCBI Gene ID8745
Uniprot	O75077
GeneID	8745;
Excitation Emission	494nm 518nm
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Background

May play a role in cell-cell and cell-matrix interactions. This is a non-catalytic metalloprotease-like protein.

Note: This product is for in vitro research use only