

## ADAM28 Antibody PE Conjugated

Catalog No: #C06055P

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	ADAM28 Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 350-400 775 derived from human ADAM28
Conjugates	PE
Target Name	ADAM28
Other Names	MDCL; MDC-L; eMDCII; ADAM 28; eMDC II; Disintegrin and metalloproteinase domain-containing protein 28; Epididymal metalloproteinase-like, disintegrin-like, and cysteine-rich protein II; Metalloproteinase-like, disintegrin-like, and cysteine-rich protein L; ADAM28; ADAM23
Accession No.	Swiss-Prot#Q9UKQ2NCBI Gene ID10863
Uniprot	Q9UKQ2
GeneID	10863;
Excitation Emission	480,565nm 578nm
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200 IF=1:50-200

## Background

May play a role in the adhesive and proteolytic events that occur during lymphocyte emigration or may function in ectodomain shedding of lymphocyte surface target proteins, such as FASL and CD40L. May be involved in sperm maturation.

Note: This product is for in vitro research use only