

ADORA2B Antibody PE Conjugated

Catalog No: #C06087P

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ADORA2B Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 145-178 332 derived from human ADORA2B
Conjugates	PE
Target Name	ADORA2B
Other Names	ADORA2; Adenosine receptor A2b; ADORA2B
Accession No.	Swiss-Prot#P29275NCBI Gene ID136
Uniprot	P29275
GeneID	136;
Excitation Emission	480,565nm 578nm
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

Receptor for adenosine. The activity of this receptor is mediated by G proteins which activate adenylyl cyclase.

Note: This product is for in vitro research use only