

GADD45B Antibody PE Conjugated

Catalog No: #C06090P

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GADD45B Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human GADD45B.
Conjugates	PE
Target Name	GADD45B
Other Names	MYD118; GADD45BETA; Growth arrest and DNA damage-inducible protein GADD45 beta; Myeloid differentiation primary response protein MyD118; Negative growth regulatory protein MyD118; GADD45B
Accession No.	Swiss-Prot#O75293NCBI Gene ID4616
Uniprot	O75293
GeneID	4616;
Excitation Emission	480,565nm 578nm
Cell Localization	Cytoplasm, Nucleus
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Background

Involved in the regulation of growth and apoptosis. Mediates activation of stress-responsive MTK1 MEKK4 MAPKKK.

Note: This product is for in vitro research use only