

CYPA Antibody FITC Conjugated

Catalog No: #C06094F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CYPA Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 130-165 165 derived from human Cyclophilin A.
Conjugates	FITC
Target Name	CYPA
Other Names	CYPA; CYPH; HEL-S-69p; Peptidyl-prolyl cis-trans isomerase A; PPIase A; Cyclophilin A; Cyclosporin A-binding protein; Rotamase A; PPIA
Accession No.	Swiss-Prot#P62937NCBI Gene ID5478
Uniprot	P62937
GeneID	5478;
Excitation Emission	494nm 518nm
Cell Localization	Cytoplasm, Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Background

PPIases accelerate the folding of proteins. It catalyzes the cis-trans isomerization of proline imidic peptide bonds in oligopeptides.

Note: This product is for in vitro research use only