

MAS1 Antibody FITC Conjugated

Catalog No: #C06101F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MAS1 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	ICC IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 70-120 325 derived from human MAS1
Conjugates	FITC
Target Name	MAS1
Other Names	MAS; MGRA; Proto-oncogene Mas; MAS1
Accession No.	Swiss-Prot#P04201NCBI Gene ID4142
Uniprot	P04201
GeneID	4142;
Excitation Emission	494nm 518nm
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Background

Receptor for angiotensin 1-7 (By similarity). Acts specifically as a functional antagonist of AGTR1 (angiotensin-2 type 1 receptor), although it up-regulates AGTR1 receptor levels. Positive regulation of AGTR1 levels occurs through activation of the G-proteins GNA11 and GNAQ, and stimulation of the protein kinase C signaling cascade. The antagonist effect on AGTR1 function is probably due to AGTR1 being physically altered by MAS1.

Note: This product is for in vitro research use only