

NARG2 Antibody FITC Conjugated

Catalog No: #C06127F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NARG2 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 240-290 982 derived from human NARG2
Conjugates	FITC
Target Name	NARG2
Other Names	BRCC1; NARG2; Little elongation complex subunit 2; Interactor of little elongator complex ELL subunit 2; NMDA receptor-regulated protein 2; ICE2; UNQ3101 PRO10100
Accession No.	Swiss-Prot#Q659A1NCBI Gene ID79664
Uniprot	Q659A1
GeneID	79664;
Excitation Emission	494nm 518nm
Cell Localization	Nucleus
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=2ug/Test IF=1:50-200

Background

Component of the little elongation complex (LEC), a complex required to regulate small nuclear RNA (snRNA) gene transcription by RNA polymerase II and III.

Note: This product is for in vitro research use only