

CD200 Antibody FITC Conjugated

Catalog No: #C06163F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD200 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 45-95 278 derived from human CD200
Conjugates	FITC
Target Name	CD200
Other Names	MRC; MOX1; MOX2; OX-2; OX-2 membrane glycoprotein; CD200; My033
Accession No.	Swiss-Prot#P41217NCBI Gene ID4345
Uniprot	P41217
GeneID	4345;
Excitation Emission	494nm 518nm
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

Costimulates T-cell proliferation. May regulate myeloid cell activity in a variety of tissues.

Note: This product is for in vitro research use only