

PRSS8 Antibody FITC Conjugated

Catalog No: #C06187F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PRSS8 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 80-130 343 derived from human PRSS8
Conjugates	FITC
Target Name	PRSS8
Other Names	CAP1; PROSTASIN; Channel-activating protease 1; Serine protease 8; PRSS8
Accession No.	Swiss-Prot#Q16651NCBI Gene ID5652
Uniprot	Q16651
GeneID	5652;
Excitation Emission	494nm 518nm
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

Possesses a trypsin-like cleavage specificity with a preference for poly-basic substrates. Stimulates epithelial sodium channel (ENaC) activity through activating cleavage of the gamma subunits (SCNN1G).

Note: This product is for in vitro research use only