

MTSS1 Antibody FITC Conjugated

Catalog No: #C06195F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MTSS1 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human MTSS1
Conjugates	FITC
Target Name	MTSS1
Other Names	MIM; MIMA; MIMB; Metastasis suppressor protein 1; Metastasis suppressor YGL-1; Missing in metastasis protein; MTSS1; KIAA429
Accession No.	Swiss-Prot#O43312NCBI Gene ID9788
Uniprot	O43312
GeneID	9788;
Excitation Emission	494nm 518nm
Cell Localization	Cytoplasm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

May be related to cancer progression or tumor metastasis in a variety of organ sites, most likely through an interaction with the actin cytoskeleton.

Note: This product is for in vitro research use only