

Wnt8A Antibody FITC Conjugated

Catalog No: #C06235F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Wnt8A Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 200-250 351 derived from human Wnt8A
Conjugates	FITC
Target Name	Wnt8A
Other Names	WNT8D; Protein Wnt-8a; Protein Wnt-8d; WNT8A
Accession No.	Swiss-Prot#Q9H1J5NCBI Gene ID7478
Uniprot	Q9H1J5
GeneID	7478;
Excitation Emission	494nm 518nm
Cell Localization	Secreted
Concentration	1mg/ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:100-500

Background

Ligand for members of the frizzled family of seven transmembrane receptors. May play an important role in the development and differentiation of certain forebrain structures, notably the hippocampus.

Note: This product is for in vitro research use only